PRESTATION

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*PARTICULAR DECONTAMINATION
AND/OR STERILIZATION*

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
Stages of decontamination - sterilization

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PARTICULAR DECONTAMINATION*

Our gloves and sleeves are decontaminated and packed in a ISO 6 clean room. The cleaning is done through surface under aqueous phase in order to obtain a classification 200 complying with the IEST STD-CC-1246D standard. Therefore, these products can be used in class 100 isolators.

DESCRIPTION

Control at reception

- Pre-cleaning (humid cloth)
- The cleaning technic is a spray of deionized H₂O filtered at 0.2 μ which match purified water quality according pharmacopeia standards. External and Internal side are decontaminated.
- Drying under laminar flow class 100 of compressed air, de-oiled and filtered at 0.22 μ .

PACKING

Standard double packing with easy opening

- 1st level class 1000 of transparent PE-type with product label - Closure by thermosealing
- 2nd level not classed transparent PE-type - Closure by thermo sealing
- Labelling : 1st level green label, product label on 1st level bag.
- Piercan will provide a COC at each order.

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GAMMA RAY STERILIZATION*

After decontamination (cleanroom laundering), the products are irradiated via gamma radiation in order to achieve sterility.

What is irradiation?

This process consists in treating the products by exposing them to gamma rays through a source of Cobalt 60 during a limited period of time.

- The treatment is secure, reliable and reproducible.
- Radio-sterilization is an environmentally friendly process.

Why to choose Gamma ray sterilization?

- The Gamma ray is particularly penetrative even at low irradiation dosages.
- This sterilization process does not mark or damage the appearance of the gauntlets in any way
- The process is not conducted at elevated temperatures and there is no detriment to the properties of the material.
- The products are sterilised inside their package.

The process we use

The bags containing the products are packaged in special cardboards (validated) and on placed into pallets. The pallets are directly exposed to a Gamma dose of between 25 to 50 KGray.

Piercan has developed a fully traceable service including the monitoring, labeling and traceability of the process methods employed, both in-house and including third party elements.

STERILIZATION VALIDATION

according to the method Vdmax 25 and norms NF EN ISO 11137 and 11737 standards

Objective and approach

Piercan has validated all stages of the process to ensure the sterility of the products after decontamination and Gamma ray sterilization. All the sterilization process by gamma irradiation of the subcontracting company has been validated.

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PARTICULAR DECONTAMINATION - STERILIZATION

Piercan products involved

- Sleeves or sleeves mounted on connectors
- ELS gloves with or without o-rings
- Dry box gloves

Piercan materials concerned

- CSM
- CSM PUR CSM
- Neoprene
- Natural rubber chlorinated

Minimum quantity per purchase order

- 30 units for gloves for dry box gloves or 30 sleeves
- 50 units ELS gloves for isolator

Upon demand

- Analysis results can be provided before and after decontamination and after sterilization (after quotation).

PIERCAN

Siège social

11 et 11bis rue Charbonnel
75013 PARIS
Tel +33 (0)1 45 88 66 27
Fax +33 (0)1 45 80 98 30
piercan@piercan.fr

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PIERCAN

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ZI Huppain - Impasse des Macareux
14520 PORT EN BESSIN
Tel +33 (0)2 31 21 73 80
Fax +33 (0)2 31 21 40 23
piercan@piercan.fr

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PIERCAN USA INC.

Filiale

180 Bosstick Blvd
San Marcos - CA 92069
Tel +1 760 599 4543
customerservice@piercanusa.com