## NOTICE D'INSTRUCTION - GANT DE BOITE A GANT EN BUTYL HAUTES PROPRIETES/EPDM CHIMIE (BHP/Eb Chim) 6/10mm MONTE SUR BAGUE SUPPORT

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## INSTRUCTIONS - GLOVEBOX GLOVES MADE OF BUTYL HIGH PROPERTIES/EPDM CHEMISTRY (BHP/Eb Chim) 6/10mm MOUNTED ON SUPPORT RING

➤ **Scope/Characteristics of the material and expiry**

- This glove can be used in nuclear industries and miscellaneous industries.
- This glove is manufactured using high-performance isobutylene polyisoprene (BHP), covered by black EPDM (Ethylene-Propylene-Diene-Monomer (Eb))
- This glove protects against certain chemical and mechanical risks, against radioactive contamination and provides protection against micro-organisms.
- It can be used in an atmosphere containing ozone. This glove does not protect against ionising radiation.
- The shelf life of this glove in its original packaging stored under the conditions described below is 3 years.

➤ **Legislation, Risk analysis and recommendations for use**

- This glove conforms to the regulation 2016/425.
- Please note, according to directive 89/656/EEC, the necessity to carry out a risk analysis relating to wearing the glove (e.g. breaking of leak tightness)
- Before each use, it is the responsibility of the user to check the quality and the integrity of the glove. The glove should be discarded if it has tears, holes or changes in surface appearance or colour that can show an alteration due to chemical products.
- For the assembly and disassembly of the gloves on the gloveboxes, please follow the instructions drafted by the safety manager and the glove port manufacturer.
- In case of a contamination, follow the instructions drafted by the safety manager.
- The gloves should be used at ambient temperature. Please contact the manufacturer for use under other temperature conditions.
- The gloves must be worn on clean and dry hands, with the nails cut short. Avoid wearing jewellery.
- It is recommended to not wear gloves when there is a risk of being caught by moving machine parts.
- During the use of unspecified chemical products (those not mentioned in the list of chemical products), please contact the manufacturer for more information.
- This glove can only be mounted if the client has a J2L system that allows receiving the glove (System not validated by IFTH during the EC type examination).

➤ **Components:**

- This glove does not contain substances in proportions in which they are known or suspected to have harmful effects on the user's health or hygiene under the foreseeable conditions of use.

➤ **Properties of the glove**

- Protection against chemical products evaluated as per **EN ISO 374-1+ A1 :2018**
- Resistance to deterioration evaluated as per **EN374-4 :2013**

Tested products		Level of permeation (palm and cuff)	Average deterioration –palm (%)	Average deterioration – cuff (%)	Pictogram
K	Sodium hydroxide 40%	6 of 6	14.1	9.4	 KMP
M	Nitric acid 65%	6 of 6	23.2	23.9	
P	Hydrogen peroxide 30%	6 of 6	6.1	18.3	

*This information does not reflect the actual duration of protection at the workplace, or the differentiation between the mixtures and pure chemical products. The chemical resistance has been evaluated under laboratory conditions using the samples collected in the palm and the cuff and only concerns the tested chemical product. It can be different if it is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the place of work can differ from the typical test conditions, depending on the temperature, abrasion and deterioration. When they are worn, the protective gloves provide less resistance to hazardous chemical products owing to the change in their physical properties. Movements, tears, friction or deterioration caused by contact with chemical products, etc. can considerably reduce the actual duration of use. For corrosive chemical products, deterioration can be the most important factor to be taken into account in the selection of chemical-resistant gloves. Before use, it is recommended to inspect the gloves in order to ensure that they have no defects or imperfections.*

- Protection against bacteria and moulds evaluated as per **EN ISO 374-5: 2016**  
Leak-tightness evaluated as per **EN374-2 :2014**

Tests	Level	Pictogram
Airtightness	Compliant	
Water-tightness	Compliant	

*The resistance to penetration has been evaluated under laboratory conditions and only pertain to the tested test specimen.*

## INSTRUCTIONS - GLOVEBOX GLOVES MADE OF BUTYL HIGH PROPERTIES/EPDM CHEMISTRY (BHP/Eb Chim) 6/10mm MOUNTED ON SUPPORT RING

- Protection against mechanical risks evaluated as per **EN388: 2016**

Tests	Level	Pictogram
Abrasion	2 of 4	 <b>2XX1X</b>
Cut by slicing	X of 5	
Tear	X of 4	
Perforation	1 of 4	
Tests Section as per EN ISO 13997 :1999	X	
X: test not carried out		

- Precision, Integrity and ozone cracking

Tests	Level
Precision <b>EN420+ A1:2009</b>	3 of 5
Integrity (pressure 30 mbar) <b>EN421-2010</b>	Compliant
Ozone cracking resistance <b>EN421-2010</b>	3 of 4

The highest figures correspond to the highest performance levels.

### ➤ Integrity EN421-2010

the “glove port diameter” and “form” lists are not exhaustive. For any requests, please contact PIERCAN.

GLOVE PORT DIAMETER (mm)	FORM	ORIGINAL LEAK-TIGHTNESS PRESSURE (mbar)
156	10	20
220	18,	20
250	85	20

### ➤ Marking

- Example of marking put on the glove and the packaging (see below)

PIERCAN LOGO


PIERCAN (1)

**CE** 0333 (2)10750 BHP/Eb Chim 8.5 E6 GJABS [XXXXX]

(3) (4) (5) (6) (7) (8) (9)

FR XXXXXX BATCH:XX XX (10)

EXPIRY:XX/XX/XXXX (11)

No.	DESCRIPTION
1	Manufacturer – In charge of release to the market
2	“CE” marking and No. of the body in charge of the annual inspection
3	Form
4	Length in mm
5	Material
6	Size
7	Thickness in 10th of mm
8	Glove mounted on support ring
9	Product Item Code
10	Case number and batch number (“FR” French manufacturing “US” American manufacturing)
11	Date of expiry in storage
12	“INFORMATION” pictogram
13	“Protection against radioactive contamination” pictogram EN 421 :2010
14	“Protection against mechanical risks” pictogram EN 388 :2016
15	“Protection against chemical risks” pictogram EN ISO 374-1 +A1 :2018
16	“Protection against bacteria and moulds” pictogram EN ISO 374-5 :2016

**INSTRUCTIONS - GLOVEBOX GLOVES**  
**MADE OF BUTYL HIGH PROPERTIES/EPDM CHEMISTRY (BHP/Eb Chim) 6/10mm**  
**MOUNTED ON SUPPORT RING**

➤ **Sizes available in the main ambidextrous models**

FORMS	Ø RDG (glove port diameter)	Ø BS (internal diameter of the support ring)	SIZES					
			7.5	8.5	9	9.5	10	10.5
10 (GJABS)	156	156	x	x		x		
18 (GJABS)	220	206		x		x		x

➤ **Suitable packaging for transport**

- The gloves should be transported in their original packaging.

➤ **Storage**

- Store in a dry place, away from light
- Storage temperature: 5 to 25°C.
- Do not store near electrical installations to avoid accelerated ageing

➤ **Cleaning and maintenance**

- If necessary, the gloves are washed with soapy water, rinsed with clear water, and dried under a low temperature air flow (<40°C) and powdered. The performance is not guaranteed by the IFTH after washing
- Do not use chemical products
- Do not use pointed or sharp objects like wire brushes, sandpaper or similar objects

➤ **Treatment of the glove**

- Non-recyclable glove
- It can be discarded as simple non-toxic waste if it is not soiled by a hazardous product
- If a glove is soiled, discard it through the appropriate channel

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