## NOTICE D'INSTRUCTION - GANT DE BOITE A GANT EN POLYURETHANE (PUR) 4/10mm MONTE SUR BAGUE SUPPORT

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## INSTRUCTIONS - GLOVEBOX GLOVES MADE OF POLYURETHANE (PUR) 4/10mm MOUNTED ON SUPPORT RING

➤ **Scope/Characteristics of the material and expiry**

- This glove can be used in the domain of life sciences (medical and pharmaceutical industries) and in nuclear industries.
- This glove is manufactured using Polyurethane.
- This glove protects against certain mechanical risks and against radioactive contamination.
- It can be used in an atmosphere containing ozone.
- This glove does not protect against ionising radiation.
- The shelf life of this glove in its original packaging stored under the conditions described below is 4 years.

➤ **Legislation, Risk analysis and recommendations for use**

- This glove conforms to the regulation 2016/425.
- Please note, according to directive 89/656/EEC, the necessity to carry out a risk analysis relating to wearing the glove (e.g. breaking of leak tightness)
- Before each use, it is the responsibility of the user to check the quality and the integrity of the glove. The glove should be discarded if it has tears, holes or changes in surface appearance or colour that can show an alteration due to chemical products.
- For the assembly and disassembly of the gloves on the gloveboxes, please follow the instructions drafted by the safety manager and the glove part manufacturer.
- In case of a contamination, follow the instructions drafted by the safety manager.
- The gloves should be used at ambient temperature. Please contact the manufacturer for use under other temperature conditions.
- The gloves must be worn on clean and dry hands, with the nails cut short. Avoid wearing jewellery.
- It is recommended to not wear gloves when there is a risk of being caught by moving machine parts.
- This glove can only be mounted if the client has a J2L system that allows receiving the glove (System not validated by IFTH during the EC type examination).

➤ **Components:**

- This glove does not contain substances in proportions in which they are known or suspected to have harmful effects on the user's health or hygiene under the foreseeable conditions of use.

➤ **Properties of the glove**

- Protection against mechanical risks evaluated as per **EN388: 2016**

Tests	Level	Pictogram
Abrasion	3 of 4	 3X22X
Cut by slicing	X of 5	
Tear	2 of 4	
Perforation	2 of 4	
Tests Section as per EN ISO 13997 :1999	X	
X: test not carried out		

- Precision, Integrity and ozone cracking

Tests	Level
Precision <b>EN420+ A1:2009</b>	5 of 5
Integrity (pressure 30 mbar) <b>EN421-2010</b>	Compliant
Ozone cracking resistance <b>EN421-2010</b>	4 of 4

The highest figures correspond to the highest performance levels.

## INSTRUCTIONS - GLOVEBOX GLOVES MADE OF POLYURETHANE (PUR) 4/10mm MOUNTED ON SUPPORT RING

### ➤ Marking

- Example of marking put on the glove and the packaging (see below)

PIERCAN LOGO

PIERCAN (1)

 0333 (2) 64750 PUR 10 E4 GJABS [XXXXX]

(3) (4) (5) (6) (7) (8) (9)

FR XXXXXX BATCH:XX XX (10)

EXPIRY:XX/XX/XXXX (11)

(12)

(13)

(14)

No.	DESCRIPTION
1	Manufacturer – In charge of release to the market
2	“CE” marking and No. of the body in charge of the annual inspection
3	Form
4	Length in mm
5	Material
6	Size
7	Thickness in 10th of mm
8	Glove mounted on support ring
9	Product Item Code
10	Case number and batch number (“FR” French manufacturing “US” American manufacturing)
11	Date of expiry in storage
12	“INFORMATION” pictogram
13	“Protection against radioactive contamination” pictogram EN 421 :2010
14	“Protection against mechanical risks” pictogram EN 388 :2016

### ➤ Sizes available in the main ambidextrous models

FORMS	Ø RDG (glove port diameter)	Ø BS (internal diameter of the support ring)	SIZES				
			8.5	9	9.5	10	10.5
54 and 64 (GJA+BS)	136	131			x	x	
10 (GJA + BS)	156	156	x		x	x	
18 (GJA + BS)	220	206	x		x		x

### ➤ Suitable packaging for transport

- The gloves should be transported in their original packaging.

### ➤ Storage

- Store in a dry place, away from light
- Storage temperature: 5 to 25°C.
- Do not store near electrical installations to avoid accelerated ageing

### ➤ Cleaning and maintenance

- If necessary, the gloves are washed with soapy water, rinsed with clear water, and dried under a low temperature air flow (<40°C) and powdered. The performance is not guaranteed by the IFTH after washing
- Do not use chemical products
- Do not use pointed or sharp objects like wire brushes, sandpaper or similar objects

### ➤ Treatment of the glove

- Non-recyclable glove
- It can be discarded as simple non-toxic waste if it is not soiled by a hazardous product
- If a glove is soiled, discard it through the appropriate channel

Attestations d'Examen UE de Type (AET) (EU type examination certifications) delivered by: I.F.T.H. (no. 0072) Avenue Guy de Collongue – 69134 ECULLY Cedex.

Monitoring organisation: AFNOR CERTIFICATION (no. 0333):11, rue Francis Pressensé-93571 LA PLAINE ST DENIS Cedex

PIERCAN: Impasse des Macareux ZI Huppain / 14520 PORT EN BESSIN FRANCE Tel.: 33 (0)2 31 21 73 80 Fax: 33 (0) 2 31 21 40 23 Email: piercan @ piercan.fr Website: www. Piercan.fr

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